



# J200

## CPAP NASAL PILLOW MASK

### TECHNICAL SPECIFICATIONS

#### Weight J200 Nasal Pillow Mask

80g incl. Headgear

#### Material

Mask Pillow                      Silicone  
Headgear                      Nylon and Elastane

#### Resistance

1.2 cmH<sub>2</sub>O @ 50 L/min.

5.0 cmH<sub>2</sub>O @ 100 L/min.

#### Dead space volume:

Size S: 70ml

Size M: 73ml

Size L: 76ml

The J200 nasal pillow mask offers patients an alternative solution to today's conventional CPAP interface. It's small, none intrusive design minimizes facial contact for optimum comfort, also offering the patient a larger field of view; ideal for patients wearing spectacles during therapy. The J200 nasal pillow mask includes 3 pillow cushion sizes ensuring the perfect fit, and a flexible light weight tube which allows over-head tube positioning.

#### FEATURES & BENEFITS:

- Discreet and comfortable by design
- 3 different pillow sizes
- 360° rotating elbow for preferred sleeping position
- Simple ergonomic 'One' piece design
- Reduced leaking, eliminating skin and eye irritation
- Comfortable headgear including over-head tube positioning
- User friendly, easy to assemble/disassemble for daily cleaning

A perfectly fitting mask is the key to successful CPAP therapy. The J200 nasal pillow mask; designed with compliance in mind!



Headgear



Nasal Pillow



Mask Frame

#### Drive DeVilbiss International

DeVilbiss Healthcare GmbH  
Kamenzer Straße 3  
68309 Mannheim, Germany

UK +44 (0) 121 / 521 31 40  
DE +49 (0) 621 / 1 78 98 0  
kontakt@devilbisshc.com

[www.drivedevilbiss-int.com](http://www.drivedevilbiss-int.com)

